


Lab Slides of Davidson® Dyes in Use



Caution: Dyes are permanent. Never apply these dyes to a living patient. Only use on excised tissue.

Dear Colleague,

Bradley Products, Inc. introduced **The Davidson Marking System®** in 1984. These high quality marking dyes continue to enjoy wide acceptance in hospitals, clinics and research facilities around the world.

These dyes are visible, reliable and safe. Their major application continues to be marking peripheral tumor or surgical margins. By applying one dye to the superior border, another to the inferior border, etc., one can identify the location of tumor or other pathology within a specimen. The dyes are also used to process multiple specimens in a single cassette. Specimens are marked with different colors, processed together, and then are all sectioned and stained on a single slide. This is both efficient and economical. Specimens retain colors and are easily distinguished.

In response to your requests, we are pleased to announce **NEW** products — **NEW Plastic Trays**, **NEW Sizes**, **NEW Sets**.

In addition, **The Davidson® Cryocup® System** has revolutionized frozen section processing and **The Davidson® Small Specimen Kit** solves the problem of tracking small specimens during histologic preparation.

Visit our web site at www.bradleyproducts.com for more information on all of the Davidson® products, including **NEW instructional videos**. As always, I appreciate your comments describing issues, uses and concerns regarding our products, and strongly encourage that communication.

Terence Davidson, M.D. — Professor of Head and Neck Surgery, UCSD

To Order:



P: 1300 778 068

E: sales@hurstscientific.com.au

www.hurstscientific.com.au



The original marking system for a variety of laboratory applications requiring orientation of tissue specimens.



Visit our **NEW Website!**

New Instructional Videos
New Shopping Cart
New Design

www.bradleyproducts.com

Product Information
For current pricing, see price sheet or visit our web site.



Bradley Products, Inc.

P.O. Box 201405 Bloomington, MN 55420
phone: 952-881-1430 toll-free: 800-325-7785 fax: 952-881-1873
email: dms@bradleyproducts.com internet: www.bradleyproducts.com

US Federal ID
41-0979462

CE IVD 00527

082211



Bradley Products, Inc.

Since 1984, supplying marking dyes to the world's medical laboratories.

The Davidson Marking SYSTEM[®]

6 Complete Sets
7 Colors 4 Sizes
Wood and Plastic Tray options

Meeting all your
margin marking needs.

Introducing **NEW** Plastic Holding Trays!



- Easy to clean, more hygienic
- Recyclable high density polyethylene plastic
- Long lasting
- Available in ALL sizes and sets
- Same price as wood trays
- 100% of evaluators switched to plastic trays after testing

Davidson Marking System[®] Complete Sets

Each set of dyes comes complete with a plastic OR wood holding tray, 50 applicator sticks and instructions.



Dye Sets including 2 oz. Davidson[®] Dyes

Item #2407-P	7 color set / plastic tray	(1 ea. green, yellow, black, red, blue, orange, violet)
Item #2407	7 color set / wood tray	(1 ea. green, yellow, black, red, blue, orange, violet)
Item #2406-P	6 color set / plastic tray	(1 ea. green, yellow, black, red, blue, orange)
Item #2406	6 color set / wood tray	(1 ea. green, yellow, black, red, blue, orange)
Item #2401-P	5 color set / plastic tray	(1 ea. green, yellow, black, red, blue)
Item #2401	5 color set / wood tray	(1 ea. green, yellow, black, red, blue)
Item #2403-P	3 color set / plastic tray	(3 ea. dyes — choose colors)
Item #2403	3 color set / wood tray	(3 ea. dyes — choose colors)

Dye Sets including 1 oz. Davidson[®] Dyes

Item #1207-P	7 color set / plastic tray	(1 ea. green, yellow, black, red, blue, orange, violet)
Item #1207	7 color set / wood tray	(1 ea. green, yellow, black, red, blue, orange, violet)

Dye Sets including 3cc Mini Davidson[®] Dyes

Item #1007-P	7 color set / plastic tray	(3 ea. green, yellow, black, red, blue, orange, violet)
Item #1007	7 color set / wood tray	(3 ea. green, yellow, black, red, blue, orange, violet)

Davidson Marking System[®] Individual Dyes

(x) below =
1 green, 2 yellow, 3 black,
4 red, 5 blue, 6 orange, 7 violet
(Example: 1163-1 is 2 oz. green dye)

Item 3408-(x)	1 ea. 8 oz. dye
Item 1163-(x)	1 ea. 2 oz. dye
Item 1101-(x)	1 ea. 1 oz. dye
Item 1013-(x)	1 jar — 40 3cc Mini dyes

3cc mini dyes are shipped as a jar of 40 3cc bottles of dye.

1 oz., 2 oz. and 8 oz. dyes are shipped as one bottle of dye.



Wood and Plastic Holding Trays

2 oz. Holding Trays

Item #1877-P	7 hole plastic tray
Item #1877	7 hole wood tray
Item #1866-P	6 hole plastic tray
Item #1866	6 hole wood tray
Item #1873-P	5 hole plastic tray
Item #1873	5 hole wood tray
Item #1833-P	3 hole plastic tray
Item #1833	3 hole wood tray

1 oz. Holding Trays

Item #1817-P	7 hole plastic tray
Item #1817	7 hole wood tray

3cc Holding Trays

Item #1807-P	21 hole plastic tray
Item #1807	21 hole wood tray

Applicator Sticks / Brush Caps



Item #8553 Applicator Sticks
5.75 inches in length. Pkg. of 500

Item #8554 Brush Caps
Brush caps fit our 2 oz. bottles.
Pkg. of 10



Cotton Tipped Applicator Sticks
6 inches in length. Pkg. of 250.

Item #8556 250 rounded end cotton tipped sticks
Item #8555 250 tapered end cotton tipped sticks

Wide Mouth Jars / Brushes



Brushes

(for use with 16 oz wide mouth jars)
Item #8557 1 inch brushes — pkg. of 25
Item #8558 2 inch brushes — pkg. of 25

Jar (Empty) and Brushes
(for marking large specimens)
Item #7399 1 (one) 16 oz. wide mouth jar
and 10 (ten) 1-inch brushes

Davidson[®] Cryocup System

Item #5364
Complete System
(4 Cryocup's and 1 Holder)

Item #5361
Cryocup[®] Holder
(For immersing in liquid nitrogen. 16.5 in. tall)

Item #5360
Cryocup[®]
(Inside diameter:
3-1/2 cm - 1-12/32 in)



Small Specimen Kit



Item #6601
Davidson[®] Small Specimen Kit
1 oz. green, red & blue,
250 cassette papers.
(not for margin marking)

Item #6605 500 Cassette Papers