

Scavenger

(Formaldehyde Neutraliser)



If you have formaldehyde or glutaraldehyde on your premises, Scavenger is a **MUST HAVE** solution for safety and comfort of staff, patients, students and visitors.

Scavenger eliminates from the air by a direct chemical reaction.

When neutralizing formaldehyde, **Scavenger** should be directed onto the spill in an easy excess quantity and once the source of vapour is neutralized the atmosphere levels should quickly drop away.

Scavenger is available in a 200ml bottle with a mist spray.

According to OH&S regulations, you should establish a procedure for leaks, spills or uncontrolled releases of hazardous substances.

Further standards state that in order to deal with leaks and spills, clean up equipment, chemicals for neutralizing or decontaminating spills and absorbent materials should be maintained at all premises on which toxic substances are kept or handled.

Workcover has suggested that it is necessary to limit exposure to formaldehyde and glutaraldehyde vapours as they are regarded as a sensitising agent as well as a possible carcinogen. In addition accidental vapour exposures can occur from various sources, with spillage's probably being the main sources. A spillage of formaldehyde can be neutralized with our **Scavenger** neutralizer.

We also have in our range the following safety spill products.

Neutralex Pouches

Neutra Wipes

Neutra-Pads

Netra-Form Spill Control Agent



Visit our website for more information and prices.

www.hurstscientific.com.au